



BGS SCIENCE ACADEMY & RESEARCH CENTER
Agalagurki, Chikkaballapura



VI Semester B.Sc.,
Paper-VII
Inorganic Chemistry Laboratory Manual



CONTENTS

1. Preparation of cuprammonium sulphate and determination of λ_{max} and hence CFSE.*
2. Gravimetric Estimation of Barium as Barium sulphate*.
3. Preparation and analysis of ferrous oxalate and its analysis*
4. Estimation of Zinc using EDTA.
5. Estimation of Nickel using EDTA and standard zinc sulphate.
6. Determination of Total Hardness of Water using EDTA.
7. Estimation of Copper in Brass.
8. Estimation of Percentage of Iron in Haematite ore using barium diphenyl amine sulphonate as an Internal indicator.
9. Determination of Calcium in the given sample of Lime stone.
10. Gravimetric Estimation of Nickel as Nickel dimethyl glyoximate.
11. Preparation of Sodium tri-oxalato ferrate(III) and estimation of iron.

*** Procedure Writing Experiments**

EXPERIMENT NO.1

AIM: Preparation of Cuprammonium sulphate complex, the determination of λ_{\max} and Crystal field splitting energy of the complex

Principle: A solution of copper sulphate reacts with excess of ammonia to form a deep blue complex Cuprammonium sulphate. The absorbance of light of this blue colored complex is measured at different wavelengths and a plot of absorbance vs wavelength, λ_{\max} is determined. From the Planck's equation $\Delta E = hc/\lambda$ crystal field stabilization energy can be calculated.

Procedure: Standard solution of copper sulphate containing $1\text{mg}/\text{cm}^3 \text{Cu}^{2+}$ is prepared by dissolving 400mg of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ crystals in a 100ml standard flask, 5cm^3 liquor ammonia is added to it and made up to the mark by adding distilled water. An intense blue colored solution of cuprammonium sulphate results. The absorbance of this solution is measured using different filters(wavelengths) in a spectrophotometer.

A plot of Absorbance Vs Wavelength is drawn, the wavelength of the corresponding the maximum absorb-ance is recorded from the graph which gives the value of λ_{\max} .

From the value λ_{\max} the Crystal field stabilization energy is calculated using the equation: $\text{CFSE} = hv = hc/\lambda_{\max}$, later CFSE for one mole can be calculated by multiplying the value using Avagadro number.

OBSERVATIONS: Reaction: $\text{Cu} + 4\text{NH}_3 \rightarrow [\text{Cu}(\text{NH}_3)_4]\text{SO}_4$

1. Mass of empty beaker = m_1 g.
2. Mass of empty beaker + $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ crystals..... = m_2 g.
3. Mass of empty $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ crystals transferred = $(m_1 - m_2)$ g.


249.5g of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ contains----- 63.5g of Copper
 $(m_1 - m_2)$ g. of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ contains ----- $63.5/249.5 \times (m_1 - m_2)$ g of Copper= 'X'g.

100 cm^3 of copper sulphate solution contains----- 'X'g. of copper
Hence 1 cm^3 of copper sulphate solution contains----- 0.01X g. of copper.

Calculation of λ_{max} :

Trial No.	λ (nm)	Absorbance or O.D.
1		
2		
3		
.		
.		
.		

A graph of O.D or Absorbance Vs λ is plotted, the peak value of wave length gives λ_{max}



From the graph: λ_{max} = maximum absorbance corresponding wavelength.

Crystal field splitting energy(CFSE): $= h\nu = hc/\lambda_{\text{max}}$; where $h = 6.634 \times 10^{-34}$ J.s,
 $c = 3 \times 10^8$ m/s.

$$\text{CFSE} = 6.634 \times 10^{-34} \times 3 \times 10^8 / \lambda_{\text{max}} = \text{'Y' Joules.}$$

CFSE for one mole of copper ions in the complex: $[\text{Cu}(\text{NH}_3)_4]_{2+}$
 $= 6.023 \times 10^{23} \times \text{'Y'} = \text{'Z' J/mol.}$

Result: λ_{max} = maximum absorbance corresponding wavelength =nm,
Crystal field splitting energy (CFSE) of the complex = 'Z' J/mol.

EXPERIMENT NO.2


AIM: TO ESTIMATE THE AMOUNT OF BARIUM PRESENT IN 250 Cm³ OF THE GIVEN SOLUTION GRAVIMETRICALLY AS BARIUM SULPHATE

PRINCIPLE: Barium is precipitated as Barium sulphate using barium chloride solution. A hot solution containing sulphate in faintly acidic medium is treated with warm 5% Barium chloride solution till the precipitation is complete. The precipitate is filtered using Whatmann No.42 filter paper, dried, ignited and weighed as Barium sulphate.

PROCEDURE: The given sulphate solution in the 250 cm³ standard flask is made up to the mark using distilled water and shaken well. 25 cm³ of the made up solution is pipetted out into a clean beaker. 1 cm³ of Con.HCl is added to it and diluted to 30 cm³ by adding 2 test tubes of distilled water. The contents are heated to boiling, 10 cm³ of hot 0.5M BaCl₂ solution is slowly added with constant stirring. The white solid is digested in hot water bath till the precipitate settles down. The completion of the precipitate is tested by adding a few drops of barium chloride solution. The precipitate is digested for half an hour. The precipitate is filtered in whatmann No.41/42 ashless filter paper. The precipitate is washed with hot water till the solid is freed from chloride ions that is the filtrate fails to give white precipitate with silver nitrate solution.

The completely dried filter paper along the sulphate precipitate is placed in a previously weighed porcelain crucible. The crucible is gently heated at first so that the paper chars and volatile matter is expelled. After charring, the temperature is raised to dull redness and burn off any carbon present. The crucible is heated to red hot for 10 to 15 minutes. The crucible is cooled in air and transferred in the desiccator. The crucible along with the precipitate is weighed accurately after 15 minutes till constant weight is obtained, its mass is recorded.

OBSERVATIONS:



1. Mass of empty crucible..... W_1 g =g
2. Mass of empty crucible + BaSO_4 precipitate = W_2 g =.....g
3. Mass of BaSO_4 precipitate = (W_2-W_1) g = (.....) =g

233g of Barium sulphate contains 137 g of Barium

W g of Barium sulphate contains $\frac{137 \times W}{233}$ g of Barium

233

= = x g of Barium

25 cm^3 of the solution contains x g of Barium =.....g

250 cm^3 of the solution contains $x \times 10$ g of Barium =g


Result: The mass of Barium present in the given solution=g

EXPERIMENT NO.3:


PREPARATION AND ANALYSIS OF FERROUS OXALATE AND ITS ANALYSIS (BOTH IRON AND OXALATE)

AIM: TO PREPARE FERROUS OXALATE AND TO ESTIMATE THE AMOUNT OF OXALATE (C₂O₄²⁻) & Fe²⁺ IN FERROUS OXALATE.

Principle: Ferrous oxalate is prepared by heating oxalic acid solution with ferrous ammonium sulphate solution acidified with dilute H₂SO₄ as:



Ferrous oxalate reacts with potassium permanganate in acidic medium during titration, ferrous is reduced to ferric ions and the oxalic acid formed is titrated against potassium permanganate solution, where oxalate ions can be estimated. Ferric ions in the solution is reduced to ferrous ions using zinc bits, the ferrous ions is titrated against the same per-manganate solution, the titre values are used to estimate the amount of ferrous ions.



Procedure: (A) Preparation of ferrous oxalate: About 7.5g of ferrous ammonium sulphate were weighed accurately and transferred into a 250cm³ beaker. About 60cm³ of distilled water and 10cm³ of H₂SO₄(2N) were added and warmed. About 37.5cm³ of 10% oxalic acid solution was added with constant stirring. The mixture was heated gently for 5 minutes. Yellow precipitate of ferrous oxalate was obtained. The precipitate is filtered and dried.

B) Analysis of ferrous oxalate:

Estimation of oxalate ion and Fe²⁺ in ferrous oxalate using standardized Potassium per-manganate solution

About 0.05g of ferrous oxalate crystals prepared were weighed accurately. It was transferred quantitatively in a clean conical flask and 25cm³ of dil. H₂SO₄(2N) was added to dissolve the same. It was heated till bubbles just appeared (55°C - 60°C). The hot solution was titrated against the std. potassium permanganate solution in the burette till the solution turned pale pink, the titre value is recorded (t₁).

About 2g of Zn dust was added to the above solution and boiled for 25 minutes. The residue obtained was dissolved in dil. H₂SO₄(2N) and filtered. It was washed with dil. H₂SO₄(2N) and the washings were added to the filtrate. The filtrate was titrated against std. KMnO₄ solution in the burette till the solution turned pale pink. the titre value is recorded (t₂).

Observations and Calculations: A: Preparation of ferrous oxalate:

Mass of ferrous oxalate prepared =g

Mass of ferrous oxalate crystals taken for analysis = 0.05g

Volume of KMnO_4 consumed for $\text{C}_2\text{O}_4^{2-}$ and Fe^{2+} = $t_1 \text{ cm}^3$ = cm^3

Volume of KMnO_4 consumed for Fe^{2+} only = $t_2 \text{ cm}^3$ = cm^3

B: (i) Estimation of oxalate ion in ferrous oxalate using standardized KMnO_4 solution

Volume of KMnO_4 consumed for $\text{C}_2\text{O}_4^{2-}$ only = $t_1 - t_2$ = cm^3 (y)

Since, 1 cm^3 $0.1 \text{ N KMnO}_4 \equiv 0.0044 \text{ g C}_2\text{O}_4^{2-}$

$y \text{ cm}^3$ $0.1 \text{ N KMnO}_4 \equiv 0.0044 \times y \text{ g C}_2\text{O}_4^{2-}$

= g (a g)

0.05g of ferrous oxalate contains (a) g of $\text{C}_2\text{O}_4^{2-}$

Hence 100g of ferrous oxalate contains: $\frac{100 \times (a)}{0.05}$ g of $\text{C}_2\text{O}_4^{2-}$

0.05

= g of $\text{C}_2\text{O}_4^{2-}$

Hence % of oxalate ions in ferrous oxalate is = %

(iii) Estimation of Fe^{2+} in ferrous oxalate using standardized KMnO_4 solution

Since, 1 cm^3 $0.1 \text{ N KMnO}_4 \equiv 0.005585 \text{ g Fe}^{2+}$

$t_2 \text{ cm}^3$ $0.1 \text{ N KMnO}_4 \equiv 0.005585 \times t_2 \text{ g Fe}^{2+}$ = g (b g)

0.05g of ferrous oxalate contains (b) g of Fe^{2+}

100g of ferrous oxalate contains $\frac{100 \times (b)}{0.05}$ g of Fe^{2+}

0.05

= g

Hence % of Ferrous ions in ferrous oxalate is = %

Result:

(a) Mass of ferrous oxalate prepared = g

(b) Percentage of $\text{C}_2\text{O}_4^{2-}$ in $\text{Fe C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$ =

(c) Percentage of Fe^{2+} in $\text{Fe C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$ =

OBSERVATIONS: Reactions: $\text{Na}_2\text{H}_2\text{Y} \rightarrow 2\text{Na}^+ + \text{H}_2\text{Y}^{2-}$
 $\text{Zn}^{2+} + \text{H}_2\text{Y}^{2-} \rightarrow \text{ZnY}^{2-} + 2\text{H}^+$
 $\text{Na}_2\text{H}_2\text{Y}$ = Disodium salt of EDTA – $\text{C}_{10}\text{H}_{14}\text{N}_2\text{Na}_2.2\text{H}_2\text{O}$

PART-I: a) Preparation of standard Zinc sulphate solution

1. Mass of weighing bottle with $\text{ZnSO}_4.7\text{H}_2\text{O}$ crystals(W_1) =.....g
2. Mass of Empty weighing bottle with the remainings (W_2) =.....g
3. Mass of $\text{ZnSO}_4.7\text{H}_2\text{O}$ crystals transferred(W_1-W_2) =.....g
Molarity of ZnSO_4 solution(M_1) = $\frac{(W_1-W_2) \times 4}{287.5}$ =M

b) Standardization of EDTA solution:

Burette: EDTA solution
Conical flask: 25 cm³ of ZnSO_4 solution + 5 cm³ of pH-10 buffer solution
Indicator: A pinch of Eriochrome black-T
End point: Wine red to blue.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of EDTA run down(cm ³)			

Concordant Value = cm³

$V_2M_2 = V_1M_1$, Molarity of EDTA solution(M_2) = $25 \times M_1 / V_2 = 25 \times \dots\dots\dots /$
EDTA ZnSO_4 =M

PART-II: Estimation of Zinc:

Burette: Standardised EDTA solution
Conical flask: 25 cm³ of ZnSO_4 solution(given) + 5 cm³ of pH-10 buffer solution
Indicator: A pinch of Eriochrome black-T
End point: Wine red to blue.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of EDTA run down(cm ³)			

Concordant Value = cm³

$V_3M_3 = V_2M_2$ Molarity of Zn^{2+} solution(M_3) = $V_2 \times M_2 / V_3 = \dots\dots\dots \times \dots\dots\dots / 25$
 Zn^{2+} EDTA =M

Mass of Zn^{2+} ions present in 250 cm³ of the given Zinc sulphate solution = Molar mass x Molarity/4 = $65.37 \times \text{Molarity} / 4$
=g

EXPERIMENT NO. 4

AIM: TO ESTIMATE THE AMOUNT/ MASS OF ZINC PRESENT IN 250 cm³ OF THE GIVEN SOLUTION USING 0.05M EDTA SOLUTION AND ZINC SULPHATE CRYSTALS

PRINCIPLE: Zinc ions in solution form coloured complex with the indicator as well as colourless complex with EDTA. Zinc-indicator complex is less stable than Zn-EDTA complex. At the end point all the free Zinc ions are complexed with EDTA. Now EDTA removes the metal forms the less stable metal-indicator complex releasing the free indicator which has a blue colour at pH-10.

PROCEDURE:

PART I : PREPARATION OF STANDARD ZINC SULPHATE SOLUTION AND STANDARDIZATION OF EDTA SOLUTION

About 3.6 g of pure zinc sulphate crystals are weighed accurately dissolved in about one test tube of dil.sulphuric acid and transferred into 250 cm³ standard flask and made up to the mark using distilled water and shaken well for uniform concentration.

25cm³ of the made up solution is pipetted out into a clean conical flask, 5 cm³ of ammonia-ammonium chloride buffer of pH-10 are added and about 2 crystals of Erio-chrome black-T indicator are added. The wine red solution is titrated against the given EDTA solution till the colour changes from wine red to blue colour. The titration is repeated for concordant values.

PART-II: ESTIMATION OF ZINC

The given Zinc sulphate solution in the standard flask is made up to the mark using distilled water and shaken well. 25cm³ of the made up solution is pipetted out into a clean conical flask, 5 cm³ of ammonia-ammonium chloride buffer of pH-10 are added and about 2 crystals of Eriochrome black-T indicator are added. The wine red solution is titrated against the given EDTA solution till the colour changes from wine red to blue colour. The titration is repeated for concordant values.

RESULT: The Amount or mass of Zinc present in 250 cm³ of the given Zinc sulphate solution is.....g.

OBSERVATIONS: Reactions: $\text{Na}_2\text{H}_2\text{Y} \rightarrow 2\text{Na}^+ + \text{H}_2\text{Y}^{2-}$
 $\text{Ni}^{2+} + \text{H}_2\text{Y}^{2-} \rightarrow \text{NiY}^{2-} + 2\text{H}^+$
 $\text{Na}_2\text{H}_2\text{Y}$ = Disodium salt of EDTA – $\text{C}_{10}\text{H}_{14}\text{N}_2\text{Na}_2 \cdot 2\text{H}_2\text{O}$

PART-I: a) Preparation of standard Zinc sulphate solution

1. Mass of weighing bottle with $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ crystals (W_1) =g
2. Mass of Empty weighing bottle with the remainings (W_2) =g
3. Mass of $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ crystals transferred ($W_1 - W_2$) =g
Molarity of ZnSO_4 solution (M_1) = $\frac{(W_1 - W_2) \times 4}{287.5}$ =M

b) Standardization of EDTA solution:

Burette: EDTA solution
Conical flask: 25 cm³ of ZnSO_4 solution + 5 cm³ of pH-10 buffer solution
Indicator: A pinch of Eriochrome black-T
End point: Wine red to blue.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of EDTA run down (cm ³)			

Concordant value = cm³

$V_2 M_2 = V_1 M_1$, Molarity of EDTA solution (M_2) = $25 \times M_1 / V_2 = 25 \times \dots\dots\dots /$
EDTA ZnSO₄ =M

PartB: Estimation of Ni²⁺ by back titration

Burette: Standard ZnSO_4 solution
Conical flask: 35ml. EDTA solution 25 cm³ of Ni^{2+} soln + 2 cm³ of $\text{NH}_3\text{-NH}_4\text{Cl}$ buffer
Indicator: Eriochrome Black-T
End-point: Blue to wine red

Trial No.	1	2	3
Final burette reading			
Initial burette reading			
Volume of ZnSO_4 run down in cm ³			

Concordant titre value = cm³

Volume of EDTA = $V_{\text{EDTA}} = (MV)_{\text{ZnSO}_4} / M_{\text{EDTA}}$

Volume of EDTA consumed;

$V_{\text{EDTA}}(\text{con}) = 35 - V_{\text{EDTA}} = (35 - \dots)$. $M_{\text{Ni}^{2+}} \times V_{\text{Ni}^{2+}} = M_{\text{EDTA}} V_{\text{EDTA}}(\text{cons})$


$M_{\text{Ni}^{2+}} = \frac{M_{\text{EDTA}} \times V_{\text{EDTA}}(\text{cons})}{V_{\text{Ni}^{2+}}} = \frac{M_{\text{EDTA}} \times V_{\text{EDTA}}(\text{cons})}{25} = \dots\dots\dots = \dots\dots\dots \text{M}$

Mass of Ni^{2+} in 250 cm³ = $\frac{M_{\text{Ni}^{2+}} \times 58.69}{4} = \dots\dots\dots = \dots\dots\dots \text{g}$.

EXPERIMENT NO. 5 Volumetric Estimation of Ni²⁺ using EDTA

Aim: Estimation of Ni²⁺ present in 250 cm³ of the given solution by back titration.

Principle: Ethylene diamine tetracetic acid (EDTA) forms complexes with a large number of cations such as Zn²⁺, Mg²⁺, Ni²⁺ etc. The disodium salt of EDTA has two replaceable hydrogen ions and may be represented as H₂Y. It reacts with M²⁺ as follows:



The above titration can be performed at pH 10 using Erichrome Black-T indicator. Since the above reaction involves liberation of H⁺, a buffer mixture has to be used to maintain the pH-10. The buffer mixture used is NH₃ – NH₄Cl. Ni²⁺ can be estimated by back titration.

Procedure: Part A: Preparation of standard ZnSO₄ solution and standardization of EDTA

About 3.59g of ZnSO₄·7H₂O is weighed accurately and transferred into a 250 cm³ standard flask. The weighed solid is dissolved in distilled water and a drop of conc. H₂SO₄ is added, the contents of the flask are mixed well and diluted upto the mark using distilled water. The molarity of ZnSO₄ solution is calculated.


25 cm³ of this solution is pipetted into a clean conical flask. 2 cm³ of NH₃ – NH₄Cl buffer and 3 drops of Erichrome Black-T indicator are added to it, then diluted with 75 cm³ (3t.t.) of water. The resulting mixture is titrated against the EDTA solution taken in the burette till the colour of the solution changes from wine-red to blue. The titration is repeated to get concordant values.

Part B: Estimation of Ni²⁺ by back titration

The solution given in the standard flask is made upto the mark using distilled water and mixed well. 25 cm³ of this solution is pipetted into a clean conical flask, the standardised EDTA solution is added from burette in excess (about 30-35 cm³). 2 cm³ of NH₃ – NH₄Cl buffer and 3 drops of Erichrome Black-T indicator are added to it, then diluted with 75 cm³ (3t.t.) of water. The resulting mixture is titrated against the EDTA solution taken in the burette till the colour of the solution changes from blue to wine-red. The titration is Repeated to get concordant values.

Result: Mass of Ni²⁺ estimated in the 250 cm³ of the given solution is _____g.

OBSERVATIONS:



PART-I: a) Preparation of standard Zinc sulphate solution

1. Mass of weighing bottle with $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ crystals (W_1) =g
2. Mass of Empty weighing bottle with the remainings (W_2) =g
3. Mass of $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ crystals transferred ($W_1 - W_2$) =g
Molarity of ZnSO_4 solution (M_1) = $\frac{(W_1 - W_2) \times 4}{287.5}$ =M

b) Standardization of EDTA solution:

Burette: EDTA solution
Conical flask: 25 cm³ of ZnSO_4 solution + 5 cm³ of pH-10 buffer solution
Indicator: A pinch of Eriochrome black-T
End point: Wine red to blue.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of EDTA run down (cm ³)			

Concordant value = cm³

$V_2 M_2 = V_1 M_1$, Molarity of EDTA solution (M_2) = $25 \times M_1 / V_2 = 25 \times \dots / \dots$
EDTA ZnSO₄ =M

PART-II: Determination of Total hardness of water

Burette: Standardised EDTA solution
Conical flask: 50 cm³ of the given sample of water + 5 cm³ of pH-10 buffer solution
Indicator: 3-4 crystals of Eriochrome black-T; End point: Wine red to blue.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of EDTA run down (cm ³)			

Concordant Value = cm³

$V_3 M_3 = V_2 M_2$, Molarity of Ca²⁺ solution (M^3) = $V_2 \times M_2 / V_3 = \dots \times \dots / 50$
(Ca²⁺) (EDTA)

=M

Hence Mass of CaCO₃ present in 250 cm³ = $M_3 \times 100 / 4 = \dots \text{g}$

Amount of CaCO₃ in 10⁶ cm³ = $\dots \times 10^6 / 250 = \dots \text{p}$

EXPERIMENT NO.6

DETERMINATION OF TOTAL HARDNESS OF OF WATER

AIM: TO ESTIMATE THE TOTAL HARDNESS OF THE GIVEN SAMPLE OF WATER USING 0.05M EDTA SOLUTION AND ZINC SULPHATE CRYSTALS

PRINCIPLE: Hardness of water is due to the presence of calcium and magnesium ions. These ions form characteristic coloured complexes with the indicator and colourless complex with the metal ions. The EDTA –metal ion complex is more stable than the metal-indicator complex.

The metal-ion indicator complex is wine red in colour while the free indicator is blue between the pH 7-11. During titration of a sample of hard water against EDTA, the colour of the solution changes from wine red to deep blue. Total hardness of water is calculated from the titre value and expressed in parts per millions of CaCO₃.

PROCEDURE:

PART I : PREPARATION OF STANDARD ZINC SULPHATE AND STANDARDIZATION OF EDTA SOLUTION


About 3.6 g of pure zinc sulphate crystals are weighed accurately dissolved in about one test tube of dil.sulphuric acid and transferred into 250 cm³ standard flask and made up to the mark using distilled water and shaken well for uniform concentration.

25cm³ of the made up solution is pipetted out into a clean conical flask, 5 cm³ of ammonia-ammonium chloride buffer of pH-10 are added and about 2 crystals of Erio-chrome black-T indicator are added. The wine red solution is titrated against the given EDTA solution till the colour changes from wine red to blue colour. The titration is repeated for concordant values.

PART II: DETERMINATION OF TOTAL HARDNESS OF WATER

50 cm³ of the given sample of hard water is pipette out into a clean conical flask, 5 cm³ of ammonia-ammonium chloride buffer of pH-10 are added and about 2 crystals of Erio-chrome black-T indicator are added. The wine red solution is titrated against the given EDTA solution till the colour changes from wine red to blue colour. The titration is repeated for concordant values.

RESULT: The total hardness of the given sample of water isppm of CaCO₃



OBSERVATIONS: PART-I: a) Preparation of standard solution of Potassium dichromate

1. Mass of weighing bottle with crystals of Potassium dichromate (W_1)=.....g
2. Mass of Empty weighing bottle with the remainings (W_2) =.....g
3. Mass of Potassium dichromate crystals transferred (W_1-W_2) =.....g
 Normality of FAS = $\frac{(W_1-W_2) \times 4}{49}$ =N

b) Standardization of Sodium thiosulphate solution:

Burette: Sodium thiosulphate solution
 Conical flask: 10 cm³ of K₂Cr₂O₇ solution + 2 cm³ of Con.HCl + 5 cm³ of 10% KI soln.
 Indicator: Freshly prepared starch
 End point: Disappearance of blue colour.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of Na ₂ S ₂ O ₃ run down(cm ³)			

Concordant Value = cm³

$V_2N_2 = V_1N_1$, Normality of Na₂S₂O₃ solution(N_2) = $25 \times N_1/V_2$

Hypo (K₂Cr₂O₇) = $25 \times \dots\dots\dots / \dots\dots\dots = \dots\dots\dots\text{N}$

PART-II: ESTIMATION OF COPPER

Burette: Standardized Sodium thiosulphate solution
 Conical flask: 10 cm³ of Cu²⁺ solution + NH₄OH + CH₃COOH + 5 cm³ of 10% KI soln.
 Indicator: Freshly prepared starch
 End point: Disappearance of blue colour.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of Na ₂ S ₂ O ₃ run down(cm ³)			

Concordant Value = cm³

$V_3N_3 = V_2N_2$, Normality Cu²⁺ solution(N_3) = $V_2N_2/V_3 = \dots\dots\dots \times \dots\dots\dots / 10$

Cu²⁺ Hypo =N

Mass of copper in 100 cm³ of the solution = $N_3 \times 63/10 = \dots\dots\dots \times 6.3$

= g. = w g

Percentage of Copper in Brass = $w \times 100/2 = \dots\dots\dots \times 50 = \dots\dots\dots$

EXPERIMENT NO. 7

Analysis of Copper in Brass

AIM: TO ESTIMATE THE PERCENTAGE OF COPPER IN THE GIVEN SAMPLE OF BRASS ALLOY(2g OF WHICH IS DISSOLVED IN 100 Cm³ OF THE GIVEN SOLUTION) USING POTASSIUM DICHROMATE CRYSTALS AND SODIUM THIOSULPHATE SOLUTION BY IODOMETRIC METHOD

PRINCIPLE: Brass is an alloy of copper and zinc with small amounts of tin, lead and iron. Brass is dissolved in 1:1 Nitric acid to obtain a solution containing cupric ions. Cupric ions react with potassium iodide in the presence of acetic acid to liberate iodine quantitatively, the liberated iodine is titrated against the standard solution of sodium thiosulphate.

PROCEDURE:

PART I: PREPARATION OF STANDARD SOLUTION OF POTASSIUM DICHROMATE AND STANDARDIZATION OF SODIUM THIOSULPHATE SOLUTION

About 0.49g of potassium dichromate crystals are weighed accurately and dissolved in distilled water, the solution is diluted to 100 cm³ in a standard flask and shaken well. 10 cm³ of this solution is pipetted out into a clean conical flask. 2cm³ of con.HCl and 5 cm³ of 10% solution of potassium iodide are added. The flask is covered with a watch glass and allowed to stand for two minutes. The liberated iodine is titrated against the given sodium thiosulphate solution taken in the burette till the solution turns greenish yellow. 1 cm³ of freshly prepared starch indicator is added and the titration is continued till the blue colour just disappears. The titration is repeated for concordant values.


PART-II: ESTIMATION OF COPPER

The given brass solution in 100 cm³ is made up to the mark using distilled water, shaken well. 10 cm³ of this solution is pipetted out into a clean conical flask. Ammonium hydroxide is added dropwise till a pale blue precipitate is obtained. Dilute acetic acid is added till the precipitate dissolves and then slight excess of acid is added. 5 cm³ of 10% solution of potassium iodide is added, then the flask is covered with a watch glass and allowed to stand for two minutes. The liberated iodine is titrated against the given sodium thio-sulphate solution taken in the burette till the solution turns to yellow. 1 cm³ of freshly prepared starch indicator is added and the titration is continued till the blue colour just disappears. The titration is repeated for concordant values.

RESULT: The percentage of Copper in the given sample of Brass is

OBSERVATIONS:

Reactions:



PART-I: a) Preparation of standard ferrous ammonium sulphate solution

1. Mass of weighing bottle with Ferrous Ammonium Sulphate crystals(W_1)=.....g
 2. Mass of Empty weighing bottle with the remainings (W_2) =.....g
 3. Mass of Ferrous Ammonium Sulphate crystals transferred(W_1-W_2) =.....g
- Normality of FAS = $\frac{(W_1-W_2) \times 4}{392}$ =N

b) Standardization of Potassium dichromate solution:

Burette: Potassium dichromate solution; **Conical flask:** 25 cm³ of FAS solution + H₃PO₄(5 cm³) + H₂SO₄(5 cm³) **Indicator:** Barium Diphenyl amine sulphonate indicator, **End point:** Green to violet.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of K ₂ Cr ₂ O ₇ (cm ³)			

Concordant Value = cm³

$$V_2N_2 = V_1N_1, \quad ,$$

(K₂Cr₂O₇) FAS

$$\text{Normality of K}_2\text{Cr}_2\text{O}_7 \text{ solution}(N_2) = 25 \times N_1 / V_2 = 25 \times \text{.....} / \text{.....} = \text{.....N}$$

PART-II: Estimation of Iron:

Burette: Potassium dichromate solution; **Conical flask:** 25 cm³ of reduced haematite ore solution + H₃PO₄(5 cm³) + H₂SO₄(5 cm³) **Indicator:** Barium Diphenyl amine sulphonate **End point:** Green to violet.

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of K ₂ Cr ₂ O ₇ (cm ³)			

Concordant Value = cm³

$$\text{Normality of Ore solution containing iron}(N_3) = V_2 \times N_2 / V_3 = \text{.....} \times \text{.....} / 25 = \text{.....} \text{ N}$$

Mass of Iron present in 250 cm³ of the given Haematite ore solution = $\frac{N_3 \times \text{Gm.Eq.Wt.}}{4}$

$$= \frac{\text{.....} \times 55.85}{4} = \text{.....g}$$

Percentage of Iron in the given sample of Haematite ore = $\frac{\text{Mass of Iron} \times 100}{2.5} = \text{.....} \times 40 = \text{.....}$

EXPERIMENT NO.8: ANALYSIS HAEMATITE ORE

AIM: DETERMINATION OF THE PERCENTAGE OF IRON IN THE GIVEN SAMPLE OF HAEMATITE ORE(2.5g of the ore is dissolved in in 250 cm³ of the given solution) USING POTASSIUM DICHROMATE SOLUTION AND FERROUS AMMONIUM SULPHATE CRYSTALS

PRINCIPLE: Haematite ore mainly contains iron having the composition Fe₂O₃ dissolved in 1:1 HCl and filtered. Ferric ions in the solution is reduced to ferrous state by stannous chloride and con.HCl, the resulting ferrous ions is estimated using potassium dichromate solution in acid medium and diphenyl amine as internal indicator.

PROCEDURE:

PART-I: PREPARATION OF STANDARD FERROUS AMMONIUM SULPHATE SOLUTION AND STANDARDIZATION OF POTASSIUM DICHROMATE SOLUTION

About of 10 g of ferrous ammonium sulphate crystals are weighed in a weighing bottle and transferred into 250 cm³ standard flask and transferred into a beaker, 1 cm³ of con. sulphuric acid is added, dissolved in water and transferred into 250 cm³ standard flask and made up to the mark using distilled water.

25 cm³ of the made up solution is pipetted out into a clean conical flask, one test tube of acidmixture(Half test tube of phosphoric acid and half test tube of con.sulphuric acid mixture) are added into it, 2 drops of diphenyl amine indicator is added and titrated against the given potassium dichromate solution taken in the burette till the green solution changes into violet. The titration is repeated for concordant values.




PART-II: ESTIMATION OF PERCENTAGE OF IRON IN THE HAEMATITE ORE

The given haematite ore solution in the standard flask is made up to the mark using distilled water and shaken well. 25 cm³ of the made up solution is pipetted out into a clean conical flask, to this solution half test of Con.HCl is added. The resulting solution is boiled, stannous chloride solution is added in drops till the solution becomes colourless

(to ensure the complete reduction 2 drops are added in excess), 2 test tubes of water are added and cooled to the room temperature. Half test tube of saturated mercuric chloride is added at a stretch, silky white precipitate is obtained(if it turns grey, the solution is rejected and the reduction is repeated). To this solution, 1 test tube of acid mixture is added(Half test tube of phosphoric acid and half test tube of con.sulphuric acid mixture) are added into it, 2 drops of diphenyl amine indicator is added and titrated against the standardized potassium dichromate solution till the colour of the solution changes green to violet. The titration is repeated for concordant readings.

- Result: 1. Amount/Mass of iron present in 250 cm³ of the given solution =g
2. The percentage of iron in the given sample of Haematite ore =

OBSERVATIONS: Reactions: $\text{CaCO}_3 + 2\text{HCl} \rightarrow \text{CaCl}_2 + \text{CO}_2 + \text{H}_2\text{O}$



PART I: STANDARDISATION OF POTASSIUM PERMANGANATE SOLUTION

1. Mass of weighing bottle with Sodium oxalate crystals(W_1) =.....g
2. Mass of Empty weighing bottle with the remainings (W_2) =.....g
3. Mass of Sodium oxalate crystals transferred(W_1-W_2) =.....g
 Normality of Sodium oxalate solution = $\frac{(W_1-W_2) \times 4}{67}$ =N

Burette: Potassium permanganate solution

Conical flask: 25 cm³ of sodium oxalate solution +2 t.t of dil.H₂SO₄

Indicator: KMnO₄ as self-indicator, End point: Permanent pale pink colour

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of KMnO ₄ (cm ³)			

Concordant Value = cm³

$$V_2N_2 = V_1N_1, ,$$

(KMnO₄) Sod.ox. Normality of KMnO₄ solution(N_2) = $25 \times N_1 / V_2 = 25 \times \dots\dots\dots /$
 =N

PART-II: ESTIMATION OF CALCIUM IN LIME STONE

Burette: Potassium permanganate solution

Conical flask: 25 cm³ of lime stone solution with calcium oxalate ppt. + 2 tt. of dil.H₂SO₄

Indicator: KMnO₄ as self-indicator; End point: Permanent pale pink colour

Trial No.	1	2	3
Final Burette Reading			
Initial Burette Reading			
Vol. of KMnO ₄ cm ³			

Concordant Value = cm³

$$V_3N_3 = V_2N_2, ,$$

Ca²⁺ (KMnO₄) Normality of Ca²⁺ solution(N_3) = $V_2 \times N_2 / V_3 = \dots\dots\dots \times \dots\dots\dots / 25$
 =N

Mass/Amount of Calcium in 250 cm³ of its solution = $N_3 \times 20 / 4 = \dots\dots\dots \times 5 = \dots\dots\dots \text{g}$

... Percentage of Calcium in the given sample of lime stone = $\dots\dots\dots \times 100 / 2 = \dots\dots\dots$

EXPERIMENT NO.9**ANALYSIS LIME STONE**

AIM: TO ESTIMATE THE PERCENTAGE OF CALCIUM IN THE GIVEN SAMPLE OF LIME STONE, 2g OF WHICH IS DISSOLVED IN 250 cm³ OF THE SOLUTION USING 0.1N POTASSIUM PERMANGANATE SOLUTION & SODIUM OXALATE CRYSTALS

PRINCIPLE: Lime stone is essentially calcium carbonate with varying amounts of magnesium carbonate. Lime stone is dissolved in 1:1 HCl, calcium is precipitated as calcium oxalate by adding a hot solution of ammonium oxalate. The precipitate is collected in quantitative filter paper which is dissolved in dil.sulphuric acid and heated, the hot solution is titrated against standardized solution of potassium permanganate.

PROCEDURE: PART I: STANDARDIZATION OF POTASSIUM PERMANGANATE SOLUTION


About 1.7g of sodium oxalate crystals is accurately weighed and transferred into 250 cm³ standard flask. It is dissolved in distilled water and made up to the mark for uniform concentration. 25 cm³ of the made up solution is pipetted out into a clean conical flask, 2 test tubes of dil.sulphuric acid is added. The contents of the flask are heated nearly to boiling, the hot solution is titrated against the given potassium permanganate solution till the appearance of permanent pale pink colour.

PART II: ESTIMATION OF CALCIUM IN THE LIME STONE

The lime stone solution given in the 250 cm³ standard flask is made up to the mark with distilled water and shaken well. 25 cm³ of this solution is pipetted out into a clean beaker, to this 2 drops of methyl orange indicator is added, it turns red. To this solution ammonium hydroxide solution is added in dropwise till it turns yellow. The about 3 cm³ of acetic acid is added till the solution becomes orange red. 4 test tubes of distilled water is added and heated nearly to boiling. To the hot solution 30 cm³ of ammonium oxalate solution is added slowly with constant stirring. The solution is made slightly alkaline by adding little amount of ammonium hydroxide solution. The white precipitate of calcium oxalate formed is allowed settle down or digested for half an hour. Completion of the precipitation is tested using little amount of ammonium oxalate. The precipitate is filtered through a quantitative filter paper(whatman No.41/42), it is washed with cold distilled water till the washings are freed from chloride ions. The precipitate is then transferred along with the filter paper into the same beaker in which the precipitation was carried out. Two test tubes full of dil.sulphuric acid is added and titrated against standardized KMnO₄ solution till permanent pale pink colour is obtained. The titration is repeated for concordancy.

RESULT: The percentage of Calcium in the given sample of lime stone is

OBSERVATIONS



1. Mass of empty and dry sintered glass crucible..... W_1 g =g
2. Mass of empty sintered glass crucible and NiDMG precipitate= W_2 g =.....g
3. Mass of NiDMG precipitate = $(W_2 - W_1)$ g = (.....) =g

288.7g of NiDMG complex contains58.7 g of Nickel

$$\begin{aligned} W \text{ g of NiDMG complex contains } & \dots\dots\dots \frac{58.7 \times W}{288.7} \text{ g of Nickel} \\ & = \dots\dots\dots = x \text{ g of Nickel} \end{aligned}$$

25 cm³ of the solution contains x g of Nickel =.....g

. . . 250 cm³ of the solution contains x X 10 g of Nickel =g

EXPERIMENT NO.10

AIM: TO ESTIMATE THE AMOUNT OF NICKEL PRESENT IN 250 Cm³ OF THE GIVEN SOLUTION GRAVIMETRICALLY AS NICKEL DIMETHYL GLYOXIMATE

PRINCIPLE: Nickel is precipitated by the addition of 1% ethanolic solution of dimethyl glyoxime to hot faintly acidic solution of the nickel salt and then by adding a slight excess of ammonia solution(freed from carbonate). The precipitate is washed with cold water and then weighed as Nickel dimethyl glyoximate after drying at 110°C - 120°C.

PROCEDURE: The given nickel salt solution in the 250 cm³ standard flask is made up to the mark using distilled water and shaken well. 25 cm³ of the made up solution is pipetted out into a clean beaker. 5 cm³ of dil.HCl is added to it and diluted to 150 cm³. The contents of the beaker are heated to 80°C and a slight excess of dimethyl glyoxime(10 cm³) reagent is added with stirring. Immediately dilute ammonia solution is added dropwise till it smells. Digest the precipitate by heating on a steam bath for 30 minutes.The completion of the precipitation is tested by adding few more drops of the reagent.The precipitate is allowed to settle for half an hour. The precipitate is filtered through a previously weighed sintered glass crucible. The precipitate is washed repeatedly with cold water until it is freed from chloride ions. Finally the precipitate is washed with few ml of rectified spirit to remove the excess reagent. The precipitate is dried at 110°C for about an hour. It is cooled in a desiccator and weighed till concordant mass is obtained.

RESULT: The amount of Nickel present in 250 cm³ of the given solution isg.

Experiment No. 11:**PREPARATION OF SODIUM TRI-OXALATO FERRATE (III) AND ESTIMATION OF IRON****AIM: TO PREPARE SODIUM TRI-OXALATO FERRATE (III) AND ESTIMATE THE AMOUNT OF IRON IN 0.2G OF SODIUM TRI-OXALATO FERRATE (III).**


Principle: Ferrous oxalate reacts with sodium oxalate and oxalic acid in presence of hydrogen peroxide gives Sodium tri-oxalato ferrate (III) complex. The complex when heated with con.sulphuric acid decomposes and liberates carbon dioxide. Ferric ions formed can be reduced to ferrous ions by heating with zinc dust in presence of a dilute acid. Ferrous ions can be estimated by titrating with standard potassium permanganate solution in acid medium. Amount of iron is calculated using the data 1dm^3 of 1N KMnO_4 is equivalent to 55.5 g of iron.

Procedure: Part I: Preparation of sodium tri-oxalato ferrate (III):

About 5g of sodium oxalate is dissolved in about 30 ml of water. To this about 3.25g of ferrous oxalate is added and gently heated. To the hot solution 12.5ml of 20 Volume hydrogen peroxide is added slowly with constant stirring. The mixture is boiled and 12.5ml of 10% oxalic acid solution is added slowly with constant stirring, The hot solution is filtered and 15ml of ethanol is added. The hot solution is filtered and 15ml of ethanol is added. The solution is cooled in a dark place when green crystals of trioxalato ferrate complex separate out. The complex is filtered and washed using 50% ethanol.

Part II : Estimation of iron:

About 1g of the complex prepared is weighed (W) and transferred into a china dish, 2 to 3ml of con.sulphuric acid is added, it is covered with a watch glass and gently heated. When the effervescence stops, it is cooled and transferred into a conical flask. About 2g of zinc dust and 10ml of dil.sulphuric acid is added. It is boiled for 20 minutes, filtered and the filtrate is titrated against the standardized decinormal potassium permanganate solution until a pale pink colour is obtained. The volume of potassium permanganate (V) is recorded and the process is repeated for concordant values. Mass of iron in 1g of the complex is calculated.

Reactions:**Observations and Calculations:**

$$\text{Mass of Iron in the complex} = 55.85 \times V \times 0.1 / 1000w$$

Where w is the mass of the complex sodium trioxalato ferrate.

- Result:** (a) Mass of sodium tri-oxalato ferrate (III) = _____ g
(b) Amount of iron in $\text{Na}_3[\text{Fe}(\text{C}_2\text{O}_4)_3]$ = _____ g

